

## 475 MP



Mobile Workstation designed  
for Operating room  
and interventional suites

$\pm 15$  Degree fingertip monitor  
tilt function

Accommodates up to 24" flat  
panel monitor

Rotates 360 degrees for views  
in portrait and landscape

Effortless vertical travel  
movement



The 475MP is designed for use in operating and interventional suites for viewing PACS images.

Vertical travel rods are constructed of stainless steel with positive locking system

Floor stand constructed of heavy gauge steel with cast aluminum top and base

Wheels mounted to rear of base

Leveling screws provided to assure floor contact and stability

Options: Keyboard/Bi-lateral mouse holder

