



Kodak
DirectView
CR 975 System

Boost CR productivity
with a workflow leader

Kodak DirectView

CR 975 System

“Drop-and-go” cassette processing, remote operations panel—and more!

Maximize throughput, productivity, and patient satisfaction with the KODAK DIRECTVIEW CR 975 System. A next-generation workflow leader, this system has been developed in direct response to customer feedback and evolving technology. The CR 975 System provides “drop-and-go” cassette processing that reduces radiographer queuing and waiting at the reader. Processed images can be sent automatically to pre-selected destinations without manual intervention. Distributed workflow is achieved with a wall-mounted KODAK DIRECTVIEW Remote

Operations Panel, allowing study functions to be performed in the examination room. As a result, radiographers are free to spend more time with each patient, and more time performing examinations. Key improvements include faster throughput, easier cassette loading, expanded software capabilities, and reduced service-access requirements for greater placement flexibility. Part of Carestream Health’s broad digital portfolio, the KODAK DIRECTVIEW CR 975 System offers the image quality, reliability, support, and service you expect from Carestream Health.





Fast throughput

- ▶ Process up to 101 cassettes an hour.
- ▶ Handle 16 cassettes at one time: up to 8 queued for processing, and 8 erased and ready for new imaging studies.
- ▶ Cassette is ready to reuse in 40 seconds.*
- ▶ Review an image in 34 seconds at a KODAK DIRECTVIEW Remote Operations Panel.*
- ▶ “Drop-and go” workflow virtually eliminates radiographer queuing at the reader and allows you to forward images directly to multiple network destinations (laser imager, diagnostic workstation, or archive) without radiographer intervention.
- ▶ Print the first image on film in 65 seconds using a KODAK DRYVIEW 6800 Laser Imager.

*For 35 x 43 cm cassette. Times for other radiography cassette sizes are the same or less. Time for image at a network device is dependent on network performance.

Enhance workflow and productivity

- ▶ Complete an exam in as few as four screens depending on system configuration.
- ▶ Save time and effort with a bar-code reader for patient/cassette/exam ID.
- ▶ Perform virtually all study functions in the examination room with a wall-mounted KODAK DIRECTVIEW Remote Operations Panel—reducing time away from the patient and improving patient satisfaction.
- ▶ Reprocess images if needed during image review—reducing the cost, time, dose, and inconvenience of repeat exposures.
- ▶ Integrate CR and DR systems together with the optional KODAK DIRECTVIEW Capture Link System to streamline workflow when a study includes CR and DR images.

Intuitive, powerful user interface

- ▶ Learn and use the system quickly with an intuitive interface designed for simple, consistent operation across all KODAK DIRECTVIEW CR and DR Systems.
- ▶ Use the same interface at the unit’s 17-inch flat-panel touch screen and at the remote operations panel.
- ▶ Pre-select up to 24 unique profiles of multiple destinations.
- ▶ Customize and preset image display parameters to suit department preferences.
- ▶ Choose the order in which patient records appear on the work list screen to suit user preference—alphabetical, oldest first, newest first.
- ▶ Take your choice of 10 languages—Chinese, English, French, German, Italian, Japanese, Norwegian, Portuguese, Spanish, or Swedish.

“KODAK DIRECTVIEW EVP Image Processing Software provides optimal viewing of structures with vastly different densities, which is a tremendous advantage.”

Dr. Donald Miller, President, Radiology Group
Marion General Hospital and Marion Area
Health Center, Marion, Ohio

Maintain high diagnostic standards

- ▶ Process images with built-in state-of-the-art KODAK DIRECTVIEW Image Processing Software.
- ▶ Take image quality and diagnostic confidence to a new level with optional KODAK DIRECTVIEW EVP Plus Software that provides extended image latitude without loss of detail contrast.
- ▶ Optimize image viewing with optional black surround/masking software.
- ▶ Improve image quality through optional automatic detection and suppression of gridline artifacts.

Image preparation and productivity features

- ▶ Improve image preparation with a configurable text box.
- ▶ Accept all study images without the need to access each individual image.
- ▶ Increase image quality and reduce retakes by magnifying an image area from one-quarter to full resolution to help identify patient position and movement.

- ▶ Reduce keystrokes by preloading procedure and exam information.
- ▶ Automatically remove patient records cancelled by the HIS/RIS from the CR database.
- ▶ Use an enhanced trauma button to quickly select a custom patient record when multiple trauma patients must be imaged.
- ▶ Promote productivity by identifying areas for improvement with an optional “reject image reason” data-entry tool.
- ▶ Optional DICOM work list management service class user software lets you obtain patient information directly from the HIS/RIS, eliminating the need for the radiographer to enter this information at the CR.



- ▶ Optional IHE scheduled workflow software automatically notifies the HIS/RIS of exam status to increase productivity with consistent patient and ordering information.

Simpler image management

- ▶ Store approximately 2,000 images online for rapid recall.
- ▶ Minimize errors with a system configuration that prohibits distribution without image approval.
- ▶ Route images automatically to multiple DICOM destinations.
- ▶ Users are alerted when patient data is incomplete or image delivery is pending.
- ▶ View the status of the last printer to which the system forwarded an image.
- ▶ Avoid image loss from power failure with standard on-board, uninterruptible power supply (UPS).



Patient data protection

- ▶ Configure the system to require entry of user name and password.
- ▶ Guard against unauthorized viewing of an unattended system with an automatic security screensaver.
- ▶ Minimize threats from worms and other malware with a built-in firewall.
- ▶ Use optional security audit log software with system login to prevent unauthorized users from accessing data stored in the CR system and track access to specific images to support privacy regulation requirements.

KODAK CARESTREAM RIS, PACS, Information Management Solutions

- ▶ Provide a totally integrated, enterprise-wide radiology suite that automates the various elements of the diagnostic exam process.
- ▶ Together, KODAK CARESTREAM RIS, PACS, and Information Management Solutions enable highly efficient diagnostic reporting and clinical review.

- ▶ Facilitate easy storage, retrieval, and distribution of a patient's entire radiology record—all through a synchronized desktop.
- ▶ Share and leverage images and information across multiple departments.
- ▶ Improve productivity, enhance patient services, and reduce network and storage infrastructure costs.

Choice of printing options

- ▶ Print CR images on film with a choice of KODAK DRYVIEW Laser Imagers. They can be configured to receive CR images automatically.
- ▶ Use a KODAK DRYVIEW 6800 Laser Imager and select from five film sizes including up to three online, all at 650 dpi resolution.
- ▶ Print true-size images using any KODAK DRYVIEW Laser Imager.
- ▶ Have flexibility to format with multiple images on a single sheet of film.

Expert network service and support

- ▶ Maximize your network throughput, performance, and reliability with Carestream Health network services.
- ▶ Count on Carestream Health as your single point of contact for expert network consulting, assessment, design, implementation, support, and monitoring.
- ▶ Check with your local Carestream Health representative for network service availability.

Remote management services

- ▶ Secure remote diagnostics and service can help resolve issues promptly and reduce service calls.
- ▶ Proactive notification to Carestream Health service centers helps optimize imager uptime.
- ▶ Remote assessment enables service personnel to arrive on-site with appropriate parts, expediting service.

ACCESSORIES AND OPTIONS

The following optional accessories and software expand the capabilities of your KODAK DIRECTVIEW CR System. Each of these is offered individually so that you can choose and purchase the options that meet your needs.

KODAK DIRECTVIEW Remote Operations Panel

Extend CR 975 System reach and functionality with this wall-mounted touch-screen panel that allows radiographers to perform most system functions away from the main unit. Network up to 10 remote operations panels to a single KODAK DIRECTVIEW CR System. Or network one remote operations panel to as many as eight CR systems.



KODAK DIRECTVIEW Total Quality Tool

Perform your own objective image tests and QC measurements with the same interface used for examinations. Test at your convenience—without the cost, scheduling, or disruption of third-party testing.



KODAK DIRECTVIEW EVP Plus Software

Take image quality and diagnostic confidence to a new level with this technological breakthrough that extends image latitude without loss of detail contrast. All raw image data are used and retained during image processing and enhancement.

KODAK DIRECTVIEW Capture Link System

Link up to five KODAK DIRECTVIEW CR and/or DR Systems to improve workflow and productivity in high-volume areas. Share cassette identification, scanning, and image review functions on linked CR systems. Enable CR cassette identification and image review at the DR console.

KODAK DIRECTVIEW CR Long-Length Imaging System

Capture long-bone computed radiography images—full leg and full spine—with this easy-to-use accessory. Fully automatic and highly reliable stitching software delivers images up to 17 inches wide by 51 inches long (43 x 129 cm), with few, if any, visible seams. A portable cassette system for supine and upright procedures produces up to 14 x 33 inches (35 x 84 cm).

KODAK DIRECTVIEW Remote Patient Data Entry Software

Increase productivity, improve workflow, and enhance patient satisfaction with remote entry of patient data. Enter demographic data for exams with a PC at the radiology reception desk or other preferred location. Enable your radiographer to spend more time performing examinations.



KODAK DIRECTVIEW Mammography Feature for CR[†]

Bring the convenience and image quality of the CR 975 System to digital mammography. This feature, coupled with KODAK DIRECTVIEW CR Mammography Cassettes and EHR Screens, delivers mammography quality images while retaining all CR system features and capabilities for general radiography examinations.

[†] Not currently available in the U.S.

