



# 480 DR Mobile Tilting Cassette Holder



Adjustable arms  
accommodate multiple  
size detectors

180 Tilt Capability

360 Degree Rotation

Channels for Stationary Grid

Unprecedented Mobility in Wireless  
Digital Radiography Imaging

	480 DR
	Tilting Mobile Cassette Holder

The 480 DR is designed to go anywhere and accommodate all flat panel wireless detectors.

**Construction:**

Mobile Stand - heavy gauge steel

Body - cast aluminum

Vertical Travel rods - Stainless steel

**Measurements:**

Base - 16.5" width x 24" depth

**Brake System:**

Located on the base of rear rubber casters

