

iCRco Scoliosis Package

“Long View Cassette and Image Stitching Software”

iCRco's Image stitching software feature QPC XSCAN32 software with extensive tools to further manipulate images and enhance workflow for the end-user.

Key Features

- Long View Cassette
- Automatic Long Bone Scoli Stitching Software
- Stitch Up to 3 Images with One Click
- Advanced Image Viewing and Manipulation Tools
- Quality Processing Edge Software

CR/DR Image Stitching Software



Long View Cassette



14"x36"

iCRco's 14"x36" Long View Cassette consists of two interlocked 14"x17" cassettes. Both Cassettes with IPs are exposed at the same time. The Image Stitching software automatically merges images together and allows you to manually manipulate images further.

IPs inside the cassettes overlap 1" to ensure capturing a complete long view image. Images can be optimized using iCRco Advanced Image Processing.

Ideal for all Scoliosis Applications.



The Innovative CR Company

2580 West 237th Street Torrance, CA 90505

Tel 310.921.9559 | Fax 310.542.7236

www.icrcompany.com